




Medizintechnik Sattler GmbH
Innovations for surgery



CATALOGUE SATTLER RETRACTOR

ABDOMINAL AND
ROCHARD- (UPPER ABDOMINAL) SYSTEM

GENERAL DESCRIPTION:

The "Sattler" retractor system is a universally designed surgical retractor system. The system is modular in construction and can be adapted and expanded any time as desired. The self-retaining modular construction makes it efficient and optimally cost-effective to use. This can reduce operation costs so it is an investment that will soon pay for itself.

SURGICAL SPECIALTIES:

- General Surgery
- Gynaecology
- Urology
- Pediatric Surgery
- Vascular Surgery
- Transplantation Surgery

BENEFITS AT A GLANCE:

- 100% Made in Germany
- Special items on customer request
- Patented clamp technology (very flat, light weight, stable)
- Articulating posts for a special situs adjustment in the deep
- Plastic coated retractor blades (non-dazzling)
- Very stable system, especially for obese people
- Quick and easy to install
- Optimization of operation costs
- No extensive training
- Practice-oriented, simple, and quick configuration
- Ideal for the most surgical areas
- Good movement opportunities, maximum of space
- Modular retractor systems
- Electrical safety
- Easy cleaning
- Certified in accordance with DIN EN ISO 13485:2008, Guideline 93/42 EWG
- Compatibility between retractors

COMPONENTS ABDOMINAL-SYSTEM



01100-62

Sterile field post, standard
with 1 clamp 20/14mm



01115-44

Support arm asymmetric with 2 clamps
14/14mm



01120-14

Frame arm (couple) for retractor blades

01115-43

Support arm, asymmetric (14 x 500 mm),
articulating with 2 clamps 14/14 mm



01115-45

Support arm for abdominal with ball system and
2 clamps 14/14 mm



COMPONENTS ROCHARD-SYSTEM



01100-63
Sterile field post with 2
clamps 20/14mm



01120-12
Cross bar Rochard, collapsible (2 hinges)



01120-29
Frame arm (couple) for retractor blades,
with joint, collapsible, for Rochard



CLAMPS & ACCESSORIES



01130-05
Quick - release clamp
for holder 14/8 mm



01130-07
Quick - release clamp for
holder 14/8 mm



01130-06
Quick - release clamp
for holder 14/8 mm



03FE000273
Flexible cold light cable, Ø
4,5/4,8 x 2650 mm with post



01130-12
Universal retaining system for optics and MIC
instruments with variable diameter Ø 2-12 mm

COMPONENTS

BLADE HOLDER



B 180
Blade holder with ball
pressure closure, 180mm



B 240
Blade holder with ball
pressure closure, 240mm



B 240/G
Blade holder with ball
pressure closure, 240mm,
articulating



B 400
Blade holder with ball
pressure closure, 400mm



B 400/G
Blade holder with ball pressure
closure, 400mm, articulating



R
G
W



R
G
W

Retractor blade Fritsch

Art. No.	Width in mm	Length in mm
0101	25	25
0102	50	40
0103	60	50
0104	80	60
0105	80	80
0107	90	105
0108	20	20
0109	30	20

Retractor blade Rochard

Art. No.	Width in mm	Length in mm
1201	80	80
1202	80	60
1203	60	50
1204	80	80
1205	80	60



R
G
W



R
G
F
SF
W

Retractor Blade Body-Wall

Art. No.	Width in mm	Length in mm
1301	80	120
1302	80	160
1303	80	200
1304	100	200

Retractor Blade Deaver slotted

Art. No.	Width in mm	Length in mm
0501	75	175
0601	60	170
0602	80	160
0603	120	150



R
G
F
SF
W



R
G
F
SF
W

Retractor blade Langenbeck

Art. No.	Width in mm	Length in mm
0201	20	50
0202	20	80
0203	30	100
0204	30	150
0205	30	190
0206	40	240
0207	40	300
0208	40	340
0209	50	125
0210	50	170
0211	50	220
0212	70	120
0213	70	190
0214	80	155
0215	80	200
0216	100	155
0217	100	200
0218	10	40
0219	10	50
0220	10	60
0221	20	40
0222	20	60
0223	30	70
0224	30	90
0225	30	110
0226	40	80
0227	40	100
0228	40	120
0229	50	100
0230	50	120

Rtractor blade Deaver

Art. No.	Width in mm	Length in mm
0401	20	50
0402	20	80
0403	30	90
0404	30	150
0405	30	190
0406	40	240
0407	40	300
0408	40	340
0409	50	125
0410	50	170
0411	50	220
0412	70	120
0413	70	190
0414	80	155
0415	80	200
0416	100	155
0417	100	200
0418	10	40
0419	10	60
0420	20	40
0421	20	60

All blade holders are available in 180mm, 240mm, 300mm and 400mm. Other blade holders or retractor blade (shape or size) are available on request.

R - Plastic coating
G - Articulating post
F - Malleable blade
SF - Semi-malleable blade
W - Interchangeable



R
G
F
SF
W



R
G
F
SF
W

Retractor blade Harrington

Art. No.	Width in mm	Length in mm
0701	70	200
0702	70	150

Retractor blade for prostatectomy

Art. No.	Width in mm	Length in mm
0801	50	210
0802	70	210
0803	70	240



R
G
F
SF
W



R
G
F
SF
W

Retractor blade for bladder

Art. No.	Width in mm	Length in mm
0901	73	170
0902	73	200

Retractor blade for urethra

Art. No.	Width in mm	Length in mm
1501	40	60
1502	40	80

All blade holders are available in 180mm, 240mm, 300mm and 400mm. Other blade holders or retractor blade (shape or size) are available on request.

R - Plastic coating
G - Articulating post
F - Malleable blade
SF - Semi-malleable blade
W - Interchangeable



R

G

W

Retractor blade for pelvis

Art. No.	Width in mm	Length in mm
1601	40/60	180
1602	40/60	240
1603	40/60	310



R

G

W

Retractor blade for renal vein

Art. No.	Width in mm	Length in mm
1001	25	150
1002	25	200



R

G

W

Retractor blade Breisky

Art. No.	Width in mm	Length in mm
0320	20	150
0321	25	90
0322	25	110
0323	30	180
0324	35	90
0325	35	110
0326	35	150
0327	40	180
0328	45	60



G


F

W

Finger Fence retractor blade, malleable

Art. No.	Width in mm	Length in mm
1101/F	90	150
1102/F	120	170
1103/F	90	90

ABDOMINAL RETRACTOR GENERAL BLADE KIT



Pos.	Art. No.	Description	QTY			
			Classic	Classic-Click	Premium	Premium-Click
1	01100-62	Sterile field post, standard with 1 clamp 20/14mm	1	1	1	1
2	01115-43	Support arm, asymmetric (14x500mm), articulating with 2 clamps 14/14 mm			1	1
3	01115-44	Support arm asymmetric with 2 clamps 14/14mm	1	1		
4	01120-14	Frame arm (couple) for retractor blades	2	2	2	2
5	01130-05	Quick -release clamp for holder 14/8 mm	5	5	7	7
6	0103 *	Retractor blade Fritsch 60 x 50 mm, post 180 mm			2	2
7	0104 *	Retractor blade Fritsch 80 x 60 mm, post 180 mm	2	2	2	2
8	0105 *	Retractor blade Fritsch 80 x 80 mm, post 180 mm	2	2	2	2
9	0107 *	Retractor blade Fritsch 90 x 105 mm, post 180 mm			2	2
10	0205 *	Retractor blade Langenbeck 30 x 190 mm, post 240 mm			1	1
11	0210/G *	Retractor blade Langenbeck 50 x 170 mm, articulating post 240 mm	1	1	1	1
12	0211/G *	Retractor blade Langenbeck 50 x 220 mm, articulating post 240 mm	1	1	1	1
13	0213 *	Retractor blade Langenbeck 70 x 190 mm, post 240 mm			1	1
14	0410 *	Retractor blade Deaver 50 x 170 mm, post 240 mm			1	1
15	0413/G *	Retractor blade Deaver 70 x 190 mm, articulating post 240 mm	1	1	1	1
16	0501 *	Retractor blade Deaver slotted 75 x 175 mm, post 240 mm	1	1	1	1
17	0702 *	Retractor blade Harrington 70 x 150 mm, post 240 mm	1	1	1	1
18	1101/FG *	Finger Fence retractor blade 90 x 150 mm malleable, articulating post 240 mm	1	1	1	1
19	1302/G *	Retractor blade Body-Wall 80 x 160 mm, articulating post 240 mm	1	1	1	1
20	03FE000273	Cold light cable elastically flexible Ø 4,5/4,8x2650 mm with post			1	1
21	01130-12	Universal retaining system for scopes/MIC- instruments, Ø 2-12 mm			1	1
22	B 180	Blade holder with ball pressure closure, 180mm		4		4
23	B 240	Blade holder with ball pressure closure, 240mm				2
24	B 240/G	Blade holder with ball pressure closure, 240mm, articulating		2		2

* For click system, Art. No. with additional letter /W

UPPER ABDOMINAL RETRACTOR GENERAL BLADE KIT



Pos.	Art.No.	Description	Qty	
			Rochard	Rochard-Click
1	01100-63	Sterile field post with 2 clamps 20/14mm	2	2
2	01120-12	Cross bar Rochard, collapsible (2 hinges)	1	1
3	01120-29	Frame arm (couple) for retractor blades, with joint, collapsible, for Rochard	2	2
4	01130-05	Quick -release clamp for holder 14/8 mm	9	9
5	1202 *	Sternal/Rochard retractor blade blade 80 x 60 mm, post 400 mm	2	2
6	1204 *	Rochard retractor blade 80 x 80 mm, post 400 mm	2	2
7	0105 *	Retractor blade Fritsch 80 x 80 mm, post 180 mm	4	4
8	0107 *	Retractor blade Fritsch 90 x 105 mm, post 180 mm	2	2
9	0210/G *	Retractor blade Langenbeck 50 x 170 mm, articulating post 240 mm	2	2
10	0211/G *	Retractor blade Langenbeck 50 x 220 mm, articulating post 240 mm	2	2
11	0213/F *	Retractor blade Langenbeck 70 x 190 mm ,malleable, post 240 mm	2	2
12	0410/G *	Retractor blade Deaver 50 x 170 mm, articulating post 240 mm	1	1
13	0413 *	Retractor blade Deaver 70 x 190 mm, post 240 mm	1	1
14	0501 *	Retractor blade Deaver slotted, 75 x 175 mm, post 240 mm	1	1
15	0702 *	Retractor blade Harrington 70 x 150 mm, post 240 mm	2	2
16	1101/FG *	Finger Fence retractor blade 90 x 150 mm, malleable, articulating post 240 mm	1	1
17	1302/G *	Retractor blade Body Wall 80 x 160 mm, articulating post 240 mm	1	1
18	1602/G *	Retractor blade for pelvis 40/60 x 240 mm, articulating post 240 mm	1	1
19	03FE000273	Cold light cable elastically flexible Ø 4,5/4,8x2650 mm with post	1	1
20	01130-12	Universal retaining system for scopes/MIC- instruments, Ø 2-12 mm	1	1
21	B 180	Blade holder with ball pressure closure, 180mm		4
22	B 240/G	Blade holder with ball pressure closure, 240mm, articulating		2
23	B 400	Blade holder with ball pressure closure, 400mm		2

* For click system, Art. No. with additional letter /W

Medizintechnik Sattler GmbH
ist zertifiziert nach:
DIN EN ISO 9001
DIN EN ISO 13485
Richtlinie 93/42 EWG

Medizintechnik Sattler GmbH

Werkstraße 9 - 07426 Königsee-Rottenbach - Germany
tel: + 49 (0) 36738 / 4 22 93 - fax: +49 (0) 36738 / 4 47 25
kontakt@medizintechnik-sattler.de - www.medizintechnik-sattler.de