

Medizintechnik Sattler GmbH

Innovations for the surgery

System description:

This universal retractor system represents the optimal support for the surgical intervention in the area of Struma, carotid, infant surgery. The system is very flexible through its innovative assembly parts and distinguishes gentle patient treatment by strain relief. The modular construction affords the individual adjustment according to your special demand. For this reason the application of the system can be realized in user-friendly and cost-effective manner.

Product features

100% Made in Germany


Minimally invasive retractors

Patient –friendly strain relief



RETRACTOR SYSTEM SATTLER
STRUMA, CAROTID, INFANT SURGERY

Retractor – Thyroid/ Carotid Blade Kit (SR-TCBK-o1)



pos.	Qty.	Item description	Part #
1	1	Sterile Field Post with one clamp 20/14 mm	01100-51
2	1	Support arm asymmetric with 1 clamp 14/14 mm 1 and 2 clamps 14/8 mm	01115-16
3	1	Frame arm (couple) for retractor blades	01120-16
4	5	Excentric Clamp 8/8 mm	01130-16
5	1	Retractor blade Langenbeck 10 x 60 mm post 180 mm	0220
6	2	Retractor blade Langenbeck 10 x 40 mm post 180 mm	0218
7	2	Retractor blade Fritsch 20 x 20 mm post 240 mm	0108
8	2	Retractor blade Richardson 21 x 25 mm post 180 mm	0329
9	2	Retractor blade Richardson 10 x 24 mm post 180 mm	0333
10	1	Retractor blade Volkmann 4-prong, semiblunt post 180 mm	1404

Thyroid/Carotid Blade Kit is suited for small incisions, especially for Thyroid and Carotid procedures.